


Psmd14 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00487

Item Number: CM0005148



Introduction

High-quality rabbit monoclonal antibody targeting human PSMD14 (UniProt O00487), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying proteasome function, protein degradation, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	PSMD14 Rabbit mAb
Remarks/Alias	PAD1; POH1; RPN11; PSMD14
Species	Human
GeneID (Human)	10213
GeneID	10213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-210 of human PSMD14 (O00487)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00487
Protein Weight	35kDa
Shipping	Ice bag