


Psmd8 Rabbit Pab For Wb, If/Icc, Elisa - P48556

Item Number: CM0028246



Introduction

Rabbit polyclonal antibody to PSMD8, a 26S proteasome non-ATPase regulatory subunit. Validated in WB, IF/ICC, and ELISA applications with cross-reactivity to human, mouse, and rat. Ideal for ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PSMD8 Rabbit pAb
Remarks/Alias	S14; p31; HIP6; HYPF; Nin1p; Rpn12; HEL-S-91n; PSMD8
Species	Human
GeneID (Human)	5714
GeneID	5714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-350 of human PSMD8 (NP_002803.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48556
Protein Weight	40kDa
Shipping	Ice bag