

Psph Rabbit Monoclonal Antibody For Wb, Elisa - P78330

Item Number: CM0007164



Introduction

High-quality PSPH rabbit monoclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, rat. Targets phosphoserine phosphatase, a key enzyme in L-serine biosynthesis. Ideal for metabolism and neurological disorder studies.

[Learn More](#)

Parameter	Value
Product Name	PSPH Rabbit mAb
Remarks/Alias	PSP; PSPHD; PSPH
Species	Human
Gene ID	5723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human PSPH (NP_004568.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78330
Protein Weight	25kDa
Shipping	Ice bag