

Ptpn22 Rabbit Pab For Wb, If/Icc, Elisa - Q9Y2R2

Item Number: CM0012883



Introduction

High-quality rabbit polyclonal antibody targeting PTPN22 (LYP), a key negative regulator of TCR signaling. Ideal for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat for studying autoimmune diseases and immune signaling.

[Learn More](#)

Parameter	Value
Product Name	PTPN22 Rabbit pAb
Remarks/Alias	LYP; PEP; LYP1; LYP2; PTPN8; PTPN22.5; PTPN22.6; PTPN22
Species	Human
GeneID	26191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human PTPN22 (NP_036543.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2R2
Protein Weight	92kDa
Shipping	Ice bag