


Ptpn5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P54829

Item Number: CM0019337



Introduction

PTPN5 (STEP) rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets a key phosphatase in synaptic plasticity and neuronal survival pathways.

[Learn More](#)

Parameter	Value
Product Name	PTPN5 Rabbit pAb
Remarks/Alias	STEP; STEP61; PTPSTEP; PTPN5
Species	Human
GeneID (Human)	84867
GeneID	84867
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-307 of human PTPN5 (NP_116170.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54829
Protein Weight	64kDa
Shipping	Ice bag