


Ptprg Rabbit Polyclonal Antibody For Wb And Elisa - P23470

Item Number: CM0016035



Introduction

High-quality PTPRG rabbit polyclonal antibody validated for WB and ELISA. Raised against recombinant human PTPRG (aa 250-450). Detects human, mouse, rat PTPRG, a receptor-type tyrosine-protein phosphatase gamma with a molecular weight of 162 kDa. Ideal for phosphatase research.

[Learn More](#)

Parameter	Value
Product Name	PTPRG Rabbit pAb
Remarks/Alias	PTPG; HPTPG; RPTPG; R-PTP-GAMMA; PTPRG
Species	Human
Gene ID	5793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-450 of human PTPRG (NP_002832.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23470
Protein Weight	162kDa
Shipping	Ice bag