


Ptpro Rabbit Polyclonal Antibody For Wb And Elisa - Q16827

Item Number: CM0028521



Introduction

Rabbit polyclonal antibody targeting human PTPRO (Receptor-type tyrosine-protein phosphatase O) validated for Western blot and ELISA. Cross-reacts with mouse. Suitable for kidney disease and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	PTPRO Rabbit pAb
Remarks/Alias	NPHS6; PTPU2; GLEPP1; PTP-OC; PTP-U2; PTPROT; R-PTP-O; PTPRO
Species	Human
Gene ID (Human)	5800
Gene ID	5800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-300 of human PTPRO (NP_109592.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q16827
Protein Weight	138kDa
Shipping	Ice bag