

Puma Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bxh1

Item Number: CM0010206



Introduction

High-specificity PUMA Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets p53 up-regulated modulator of apoptosis, key for cancer and cell death research.

[Learn More](#)

Parameter	Value
Product Name	PUMA Rabbit mAb
Remarks/Alias	JFY1; PUMA; JFY-1
Species	Human
GeneID (Human)	27113
Immunogen	Synthetic Peptide (amino acids 100-193 of human PUMA, NP_055232.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXH1
Protein Weight	18kDa
Shipping	Ice bag