

Puma Rabbit Pab For Wb - Q96Pg8/Q9Bxh1

Item Number: CM0020025

Introduction

PUMA Rabbit pAb validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Targets multiple PUMA isoforms (10-26 kDa) for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PUMA Rabbit pAb
Remarks/Alias	JFY1; PUMA; JFY-1
Species	Human
Gene ID (Human)	27113
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-131 of human PUMA (NP_001120713.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PG8/Q9BXH1
Protein Weight	10kDa/14kDa/20kDa/26kDa
Shipping	Ice bag