

Puromycin Rabbit Monoclonal Antibody - Validated For Wb, Ihc, If, Ip, Flow Cytometry, And Elisa

Item Number: CM0000066

Introduction

This puromycin rabbit monoclonal antibody, raised against puromycin chemical compounds, is ideal for detecting puromycin incorporation in protein synthesis assays. Validated for WB, IHC-P, IF/ICC, IP, FC (intra), and ELISA; species-independent reactivity.

[Learn More](#)

Parameter	Value
Product Name	Puromycin Rabbit mAb
Gene ID	CAS:58-58-2
Immunogen	Chemical compounds corresponding to puromycin.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IHC-P, IF/ICC, IP, FC (intra), ELISA
Cross Reactivity	Species independent
Shipping	Ice bag