

Puromycin Rabbit Pab - Cas:58-58-2

Item Number: CM0019707

Introduction

Puromycin Rabbit pAb is a rabbit polyclonal antibody developed against puromycin. Suitable for WB and ELISA with species-independent cross reactivity, this antibody serves as a tool for studying protein synthesis and puromycin-based selection.

[Learn More](#)

Parameter	Value
Product Name	Puromycin Rabbit pAb
Species	Human
Gene ID (Human)	CAS:58-58-2
Gene ID	CAS:58-58-2
Immunogen	Chemical compounds corresponding to Puromycin
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag