


Pygb Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P11216

Item Number: CM0029046



Introduction

Rabbit polyclonal antibody against human PYGB (brain glycogen phosphorylase). Suitable for WB, IF/ICC, IP, and ELISA. Detects human, mouse, and rat PYGB with a predicted molecular weight of 97 kDa. Ideal for glycogen metabolism and brain energy homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	PYGB Rabbit pAb
Remarks/Alias	GPBB; PYGB
Species	Human
Gene ID (Human)	5834
Gene ID	5834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 544-843 of human PYGB (NP_002853.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11216
Protein Weight	97kDa
Shipping	Ice bag