


Rab10 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P61026

Item Number: CM0006944



Introduction

High-quality RAB10 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat RAB10. Ideal for intracellular trafficking and membrane transport research.

[Learn More](#)

Parameter	Specification
Product Name	RAB10 Rabbit mAb
Remarks/Alias	RAB10; ras-related protein Rab-10; [KD Validated] RAB10
Species	Human
Gene ID (Human)	10890
Gene ID	10890
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 101-200 of human RAB10 (NP_057215.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61026
Protein Weight	23kDa
Shipping	Ice bag