

Rab34 Rabbit Polyclonal Antibody For Wb, Elisa - P0Di83/Q9Bzg1

Item Number: CM0014673

Introduction

Rabbit polyclonal antibody against human RAB34, validated for WB and ELISA. Detects multiple isoforms in human, mouse, and rat samples. Suitable for nucleolar function and alternative splicing studies.

[Learn More](#)

Parameter	Value
Product Name	RAB34 Rabbit pAb
Remarks/Alias	RAH; NARR; RAB39; RAB34
Species	Human
Gene ID (Human)	83871
Gene ID	83871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-259 of human RAB34 (NP_114140.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0Di83/Q9BZG1
Protein Weight	21kDa/28kDa/29kDa
Shipping	Ice bag