


Rab3D Rabbit Pab For Wb, Elisa - O95716

Item Number: CM0017225



Introduction

RAB3D Rabbit pAb (O95716) is a polyclonal antibody validated for WB and ELISA, cross-reactive with human and mouse. The predicted band is 24kDa. Ideal for studying Rab-3D in membrane trafficking and exocytosis.

[Learn More](#)

Parameter	Value
Product Name	RAB3D Rabbit pAb
Remarks/Alias	GOV; D2-2; RAB16; RAD3D; RAB3D
Species	Human
GeneID (human)	9545
GeneID	9545
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-219 of human RAB3D (NP_004274.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95716
Protein Weight	24kDa
Shipping	Ice bag