


Rab4A Rabbit Pab For Wb And Elisa - P20338

Item Number: CM0029379



Introduction

RAB4A rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects RAB4A in human, mouse, and rat samples. Target is a small GTPase regulating intracellular membrane trafficking and protein transport. Ideal for endosomal pathway studies.

[Learn More](#)

Parameter	Value
Product Name	RAB4A Rabbit pAb
Remarks/Alias	RAB4; HRES1; HRES-1; HRES-1/RAB4; RAB4A
Species	Human
Gene ID	5867 (human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human RAB4A (NP_004569.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20338
Protein Weight	24kDa
Shipping	Ice bag