

Rabbit Anti-4-Hydroxynonenal (4-Hne) Monoclonal Antibody For Intracellular Flow Cytometry

Item Number: CM0000052

Introduction

-conjugated rabbit monoclonal antibody targeting 4-Hydroxynonenal (4-HNE), a key lipid peroxidation marker. Validated for intracellular flow cytometry with species-independent reactivity. Supplied as an IVD raw material for diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-4-Hydroxynonenal mAb
Remarks/Alias	4-HNE
Immunogen	Chemical compounds corresponding to 4-Hydroxynonenal.
Source	Rabbit
Category	Other Modified Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag