

Abflo 488 Rabbit Anti-4-Hydroxynonenal (4-Hne) Monoclonal Antibody For Intracellular Flow Cytometry

Item Number: CM0000027

Introduction

Fluorescent ABflo 488-conjugated rabbit monoclonal antibody targeting the lipid peroxidation marker 4-Hydroxynonenal (4-HNE). Ideal for intracellular flow cytometry (FC intra) with species-independent reactivity, supporting oxidative stress research.

[Learn More](#)

Parameter	Detail
Product Name	Rabbit anti-4-Hydroxynonenal mAb
Remarks/Alias	4-HNE
Immunogen	Chemical compounds corresponding to 4-Hydroxynonenal
Source	Rabbit
Category	Other Modified Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag