

Rabbit Anti-Atf5 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y2D1

Item Number: CM0008457

Introduction

Rabbit anti-ATF5 monoclonal antibody for Western blot, IHC-P, and ELISA applications. Validated in human, mouse, and rat, this antibody targets the transcription factor ATF5, critical for cell survival and differentiation.

[Learn More](#)

Parameter	Value
Product Name	ATF5 Rabbit mAb
Remarks/Alias	ATFX; HMFN0395; ATF5
Species	Human
GeneID (Human)	22809
GeneID	22809
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 183-282 of human ATF5 (Q9Y2D1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2D1
Protein Weight	31kDa
Shipping	Ice bag