

Rabbit Anti-Calca Polyclonal Antibody For Wb And Elisa Research Applications - P01258 / P06881

Item Number: CM0015934



Introduction

Rabbit polyclonal antibody targeting human CALCA (calcitonin), validated for Western blot and ELISA, cross-reacts with mouse and rat. Ideal for studying calcium metabolism and hormone signaling.

[Learn More](#)

Parameter	Value
Product Name	CALCA Rabbit pAb
Remarks/Alias	CT, KC, PCT, CGRP, CALC1, CGRP1, CGRP-I, CGRP-alpha, CALCA
Species	Human
GeneID (Human)	796
GeneID	796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-141 of human CALCA (NP_001029124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01258 / P06881
Protein Weight	13kDa / 15kDa
Shipping	Ice bag