

Rabbit Anti-Cat Cd27 Monoclonal Antibody - Cd27

Item Number: CM0010137

Introduction

-conjugated rabbit monoclonal antibody against cat CD27 for flow cytometry applications. Immunogen spans amino acids 20-191 of feline CD27; no cross-reactivity with other species reported. Ideal for feline immunological studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Cat CD27 mAb
Species	Cat
Gene ID (Human)	939
Gene ID	101091633
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 20-191 of cat CD27 (XP_011282025.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cat
Shipping	Ice bag