

Rabbit Anti-Cat Cd8 Alpha Monoclonal Antibody For Fc - P41688

Item Number: CM0004404



Introduction

-conjugated rabbit monoclonal antibody against cat CD8 alpha, validated for flow cytometry. Ideal for immunophenotyping and T-cell research in feline models. Target specificity to UniProt P41688.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Cat CD8 alpha mAb
Remarks/Alias	CD8A
Species	Cat
GeneID (Human)	925
GeneID	493799
Immunogen	Recombinant protein corresponding to amino acids 25-186 of cat CD8 alpha (NP_001009843.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cat
SWISS	P41688
Protein Weight	32-34 kDa
Shipping	Ice bag