

Rabbit Anti-Cat Cd8 Alpha Monoclonal Antibody For Flow Cytometry - P41688

Item Number: CM0004508



Introduction

Recombinant rabbit monoclonal antibody against feline CD8 alpha, suitable for flow cytometry (FC). Facilitates the study of CD8+ T-cell and NK cell responses in cat immunology research. Cross-reacts specifically with cat.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Cat CD8 alpha mAb
Remarks/Alias	CD8A
Species	Cat
GeneID (Human)	925
GeneID	493799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-186 of cat CD8 alpha (NP_001009843.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cat
SWISS	P41688
Protein Weight	32-34kDa
Shipping	Ice bag