

Rabbit Anti-Cc2D1A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P1N0

Item Number: CM0024841

Introduction

Rabbit polyclonal antibody against human CC2D1A, validated for WB, IHC-P, ELISA, and cross-reactive with mouse and rat. Suitable for research on gene regulation, synaptic function, and PI3K signaling.

[Learn More](#)

Parameter	Specification
Product Name	CC2D1A Rabbit pAb
Remarks/Alias	Lgd2; MRT3; TAPE; Aki-1; FREUD-1; Freud-1/Aki1; CC2D1A
Species	Human
GeneID (Human)	54862
GeneID	54862
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CC2D1A (NP_060191.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P1N0
Protein Weight	104kDa
Shipping	Ice bag