


Rabbit Anti-Cd14 Monoclonal Antibody For Flow Cytometry - P08571

Item Number: CM0004334



Introduction

-conjugated rabbit anti-CD14 monoclonal antibody for flow cytometry. Detects CD14 on human and cynomolgus monocytes/macrophages. LPS coreceptor; critical for innate immunity research. Molecular weight: 40 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID (Human)	929
GeneID	929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag