


Rabbit Anti-Cd62E/E-Selectin Monoclonal Antibody For Wb, Fc, Elisa - P16581

Item Number: CM0006597



Introduction

Human CD62E/E-Selectin rabbit monoclonal antibody validated for WB, FC, ELISA. Recombinant immunogen, detects 66kDa band. Suitable for studying leukocyte adhesion and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CD62E/E-Selectin Rabbit mAb
Remarks/Alias	SELE; CD62E; ELAM; ELAM1; ESEL; LECAM2; E-selectin
Species	Human
GeneID (Human)	6401
GeneID	6401
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human CD62E/E-Selectin (NP_000441.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P16581
Protein Weight	66kDa
Shipping	Ice bag