


Rabbit Anti-Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0012511



Introduction

Rabbit monoclonal antibody against human CD8a, conjugated to , for flow cytometry. Validated on human and cynomolgus monkey samples. Ideal for immune cell profiling in immunology and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8a (NP_001759.3))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag