

Rabbit Anti-Cdc42 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P60953

Item Number: CM0017460



Introduction

Rabbit polyclonal antibody against CDC42 (Cell division control protein 42 homolog). Validated in WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Targets small GTPase regulating actin dynamics, cell polarity, and migration.

[Learn More](#)

Parameter	Value
Product Name	CDC42 Rabbit pAb
Remarks/Alias	TKS; G25K; CDC42Hs; CDC42
Species	Human
GeneID (Human)	998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-191 of human CDC42 (NP_001782.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60953
Protein Weight	21kDa
Shipping	Ice bag