

Rabbit Anti-Ddddk-Tag Monoclonal Antibody - Ddddk-Tag

Item Number: CM0010690

Introduction

-conjugated rabbit monoclonal antibody targeting DDDDK-tag for intracellular flow cytometry. Species-independent tag recognition enables flexible experimental designs across various model systems.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-DDDDK-Tag mAb
Remarks / Alias	DDDDK; DDDDK tag; DDDDK-tag
Immunogen	Synthetic peptide: A synthetic peptide corresponding to DDDDK tag
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag