

Rabbit Anti-Dog Cd19 Monoclonal Antibody - Cd19

Item Number: CM0004723

Introduction

-conjugated rabbit monoclonal antibody targeting dog CD19, a B-cell surface antigen. Designed for flow cytometry (FC) applications in canine immunology research. Cross-reacts with dog. Suitable for direct immunofluorescence staining.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Dog CD19 mAb
Remarks/Alias	CD19
Species	Dog
Gene ID (Human)	930
Gene ID	607898
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
Shipping	Ice bag