

## Rabbit Anti-Dog Cd19 Monoclonal Antibody - Cd19

Item Number: CM0004626

### Introduction

-conjugated rabbit monoclonal antibody against dog CD19 for direct flow cytometry (FC) analysis of canine B cells. Recombinant immunogen, specific to dog CD19.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Dog CD19 mAb                                                                                          |
| Alias            | CD19                                                                                                              |
| Species          | Dog                                                                                                               |
| Human Gene ID    | 930                                                                                                               |
| Gene ID          | 607898                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-299 of dog CD19 (XP_013969991.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Dog                                                                                                               |
| Shipping         | Ice bag                                                                                                           |