

Rabbit Anti-Dog Cd27 Monoclonal Antibody For Flow Cytometry - Cd27

Item Number: CM0010234

Introduction

-conjugated rabbit monoclonal antibody targeting dog CD27 (TNFRSF7). Recombinant immunogen (aa 20-215) designed for flow cytometry. Suitable for canine research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Dog CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Dog
Gene ID (Human)	939
Gene ID	611674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-215 of dog CD27 (XP_038295140.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
Protein Weight	31kDa
Shipping	Ice bag