

Rabbit Anti-Dog Cd4 Monoclonal Antibody For Flow Cytometry - Canine Cd4

Item Number: CM0008268

Introduction

-conjugated rabbit monoclonal antibody targeting canine CD4, validated for flow cytometry (FC). Recombinant immunogen (aa 26-400). Ideal for dog lymphocyte phenotyping and immunological research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Dog CD4 mAb
Remarks/Alias	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Laurasiatheria; Carnivora; Caniformia; Canidae; Canis
Species	Dog
Gene ID (Human)	920
Gene ID	403931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-400 of dog CD4 (NP_001003252.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
Shipping	Ice bag