

Rabbit Anti-Dog Cd4 Monoclonal Antibody - Np_001003252.1

Item Number: CM0004453

Introduction

Rabbit anti-Dog CD4 Monoclonal Antibody for flow cytometry (FC) applications. Recombinant protein immunogen, reacts with dog CD4. Ideal for canine immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Dog CD4 mAb
Remarks/Alias	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Laurasiatheria; Carnivora; Caniformia; Canidae; Canis
Species	Dog
Gene ID (Human)	920
Gene ID	403931
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
Shipping	Ice bag