

Rabbit Anti-Dog Cd8A Monoclonal Antibody For Flow Cytometry - P33706

Item Number: CM0010984

Introduction

High-quality -conjugated rabbit monoclonal antibody against dog CD8A for flow cytometry. Targets T-cell surface glycoprotein CD8 alpha chain, essential for MHC class I interaction and cytotoxic T-cell responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Dog CD8A mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha; CD8A
Species	Dog
GeneID (Human)	925
GeneID	403157
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-139 of dog CD8A (NP_001002935.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
SWISS	P33706
Protein Weight	26 kDa
Shipping	Ice bag