

Rabbit Anti-Fitc/5-Fam/6-Fam Mab - Fitc

Item Number: CM0008805

Introduction

-conjugated rabbit monoclonal antibody targeting FITC, 5-FAM, and 6-FAM dyes. Species-independent reactivity and validated for flow cytometry, suitable for fluorescent probe detection in research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-FITC/5-FAM/6-FAM mAb
Remarks/Alias	FITC
Immunogen	Chemical compounds Chemical/ Small Molecule corresponding to FITC/5-FAM/6-FAM.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag