

Rabbit Anti-Fitc/5-Fam/6-Fam Mab - Fitc

Item Number: CM0012252

Introduction

-conjugated rabbit monoclonal antibody specifically recognizing FITC, 5-FAM, and 6-FAM small molecules. This species-independent reagent is designed for flow cytometry, supporting immunofluorescence detection in diagnostic assay development and biomedical research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-FITC/5-FAM/6-FAM mAb
Remarks/Alias	FITC
Species	Rabbit
Immunogen	Chemical compounds corresponding to FITC/5-FAM/6-FAM
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag