

Rabbit Anti-Fitc/5-Fam/6-Fam Mab - Fitc

Item Number: CM0004082

Introduction

Rabbit monoclonal anti-FITC/5-FAM/6-FAM antibody () for flow cytometry. Recognizes FITC, 5-FAM, and 6-FAM with species-independent reactivity. Useful for detecting FITC conjugates in immunofluorescence and IVD research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-FITC/5-FAM/6-FAM mAb
Remarks/Alias	FITC
Immunogen	Chemical compounds
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag