

Rabbit Anti-Gap43 Polyclonal Antibody For Wb, Ihc-P, Elisa - P17677

Item Number: CM0022235

Introduction

Rabbit polyclonal antibody recognizing GAP43 (neuromodulin), optimized for Western blot, IHC-P, and ELISA. Suitable for human, mouse, and rat neural tissue samples.

[Learn More](#)

Parameter	Value
Product Name	GAP43 Rabbit pAb
Remarks/Alias	B-50; PP46; GAP-43; GAP43
Species	Human
Gene ID (Human)	2596
Gene ID	2596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human GAP43 (NP_002036.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P17677
Protein Weight	25kDa
Shipping	Ice bag