

Rabbit Anti-Gnai3 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P08754

Item Number: CM0015154



Introduction

Rabbit polyclonal antibody against human GNAI3, validated for Western blot, IF/ICC, IP, and ELISA. Recognizes human and mouse GNAI3; suitable for studying G protein signaling and cell division.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GNAI3 Rabbit pAb
Remarks/Alias	87U6; HG1A; ARCND1; ARCODS; I3
Species	Human
GeneID	2773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-354 of human GNAI3 (NP_006487.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08754
Protein Weight	41kDa
Shipping	Ice bag