

## Rabbit Anti-Goat IgG (Fab) Fragment Mab - Goat IgG (Fab) Fragment

Item Number: CM0006154

### Introduction

The Rabbit anti-Goat IgG (Fab) fragment mAb is a rabbit monoclonal antibody specifically recognizing goat IgG Fab, suitable for flow cytometry applications in immunoassays and diagnostics.

[Learn More](#)

| Parameter        | Value                                                      |
|------------------|------------------------------------------------------------|
| Product Name     | Rabbit anti-Goat IgG (Fab) fragment mAb                    |
| Immunogen        | Recombinant protein of rabbit anti-Goat IgG (Fab) fragment |
| Source           | Rabbit                                                     |
| Category         | Monoclonal Antibodies                                      |
| Application      | FC                                                         |
| Cross Reactivity | Goat                                                       |
| Shipping         | Ice bag                                                    |