

Rabbit Anti-Goat Igg (Fab) Fragment Mab - Goat Igg (Fab) Fragment

Item Number: CM0004024

Introduction

-conjugated rabbit monoclonal antibody targeting goat IgG Fab fragment. Validated for flow cytometry applications. Recombinant immunogen ensures high specificity and lot-to-lot consistency. Suitable for immunoassays requiring anti-goat IgG detection in research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Goat IgG (Fab) fragment mAb
Species	Goat
Immunogen	Recombinant protein of Goat IgG (Fab) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag