

Rabbit Anti-Goat IgG (Fab) Fragment Monoclonal Antibody For Flow Cytometry - Goat IgG Fab Fragment

Item Number: CM0006402

Introduction

Rabbit monoclonal antibody conjugated to targeting the goat IgG Fab fragment. Designed for flow cytometry (FC) applications, offering species-specific detection of goat IgG in various immunoassays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Goat IgG (Fab) fragment mAb
Species	Goat
Immunogen	Recombinant protein of Goat IgG (Fab) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag