

Rabbit Anti-Goat IgG (Fc) Fragment Mab

Item Number: CM0007327

Introduction

-conjugated rabbit monoclonal antibody targeting goat IgG Fc fragment, validated for flow cytometry (FC). Ideal for immunoassays and secondary detection in goat sample systems.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of Goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag