

Rabbit Anti-Goat IgG (Fc) Fragment Mab

Item Number: CM0010190

Introduction

-conjugated rabbit monoclonal antibody targeting goat IgG (Fc) fragment, suitable for flow cytometry applications. Designed for specific detection of goat IgG in research and diagnostic assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit Anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag