

Rabbit Anti-Goat IgG (Fc) Fragment Mab - Goat IgG (Fc) Fragment

Item Number: CM0006533

Introduction

-conjugated rabbit monoclonal antibody specific for goat IgG Fc fragment. Validated for flow cytometry, suitable for detecting goat immunoglobulins in research and diagnostic assays.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of Goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag