

## Rabbit Anti-Goat IgG (Fc) Fragment Mab

Item Number: CM0007217

### Introduction

conjugated Rabbit monoclonal antibody targeting Goat IgG (Fc) fragment for flow cytometry applications. Specific for Goat Fc region, minimizing cross-reactivity. Ideal for multi-color immunofluorescence studies.

[Learn More](#)

| Parameter        | Value                                         |
|------------------|-----------------------------------------------|
| Product Name     | Rabbit anti-Goat IgG (Fc) fragment mAb        |
| Species          | Goat                                          |
| Immunogen        | Recombinant protein of Goat IgG (Fc) fragment |
| Source           | Rabbit                                        |
| Category         | Monoclonal Antibodies                         |
| Application      | Flow Cytometry (FC)                           |
| Cross Reactivity | Goat                                          |
| Shipping         | Ice bag                                       |