

Rabbit Anti-Goat IgG (Fc) Fragment Monoclonal Antibody

Item Number: CM0005290

Introduction

-conjugated rabbit monoclonal antibody targeting goat IgG Fc fragment for flow cytometry. Ideal for detecting goat IgG in immunological assays. Cross-reacts with goat IgG.

[Learn More](#)

Parameter	Value
Product Name	Rabbit Anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag