

Rabbit Anti-Hnf4A Polyclonal Antibody For Wb, Elisa - Hnf4A

Item Number: CM0013448

Introduction

Rabbit polyclonal antibody targeting HNF4a (HNF4alpha), validated for Western blot and ELISA. Cross-reacts with mouse. Useful for research into liver-specific gene regulation and hepatocellular function.

[Learn More](#)

Parameter	Value
Product Name	HNF4a Rabbit pAb
Remarks/Alias	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha; HNF4a
Species	Mouse
Gene ID (Human)	3172
Gene ID (Mouse)	15378
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human HNF4a (NP_000448.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
Protein Weight	51kDa/52kDa
Shipping	Ice bag