


Rabbit Anti-Hpcal1 Polyclonal Antibody For Wb And Elisa - P37235

Item Number: CM0024317



Introduction

High-quality rabbit polyclonal antibody against HPCAL1 (Hippocalcin-like protein 1), validated by KO, for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HPCAL1 Rabbit pAb
Remarks/Alias	BDR1; HLP2; VILIP-3; L1
Species	Human
GeneID (Human)	3241
GeneID	3241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human HPCAL1 (NP_002140.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37235
Protein Weight	22kDa
Shipping	Ice bag